

Storage of CRRT Blood with Pediatric Patients on ECMO During Transports to Reduce Blood Exposure

James Neal, CCP FPP
Assistant Professor of Surgery
Mayo Clinic
Rochester, MN
neal.james@mayo.edu



Disclosures

- None



Goals

- Show the need for a system to save blood on ECMO for pediatric and neonate patients
- Review simple math steps to repurpose a standard CPD bag
- Other uses of this technique in ECMO



Background

- Using CRRT in pediatric and neonate patients while on ECMO is becoming more published
- CRRT circuit can contain a large percentage of total patient blood
- Transports do happen with these patients and causes an issue of what to do with the CRRT circuit
- Excessive or unneeded transfusions are not helpful in pediatric and neonate patients



Adapting the autologous blood storage bag

- CPD bag has enough citrate to anticoagulant 450 mL of whole blood
- Pediatric CRRT circuit holds a residual volume of 165 mL
- CPD bag has 63 mL of CPD so the amount in the bag must be reduced precisely

MAIO CLINIC

Benefits of CPD as anticoagulant

- Reversible with just Ca administration
- Will not cause a long term increase patient coagulation values
- Can safely anticoagulant blood for long periods of time
- Used in many areas of the hospital including autologous blood removal during cardiac surgery

MAIO CLINIC

Methods (Shortened)

- All of the CPD is removed from the bag and then only the calculated amount is added back in
- Clip off long end of tubing that holds CPD
- Need to correctly measure blood entering bag via counter on CRRT machine or syringes if needed

MAIO CLINIC

Solving for the amount of CPG needed

No Citrate Infusing into CRRT circuit example:

$$\frac{63\text{mL}}{450\text{mL}} \times \frac{x\text{mL}}{165\text{mL}} = 23\text{mL CPD returned to bag}$$

Solve for x using a ratio

MAIO CLINIC

What if using citrate in CRRT circuit

Citrate
Infusing
into
CRRT
circuit
example:

$$\frac{20\text{mg/dL}^*}{50\text{mg/dL}^{**}} \times 23\text{mL} = 9.2\text{mL CPD returned to bag}$$

*CRRT Post-Filter Calcium Level
**Assumed Normal Calcium Level

MAVO CLINIC



MAVO CLINIC

Other Thoughts to Consider

- Patient may need small bolus amounts of CaCl when CRRT is restarted
- This will treat the hypocalcemia caused by the CRRT circuit volume
- Make sure to label this bagged product within established guidelines

MAVO CLINIC

Results and Conclusions

- At our center, storage of CRRT circuit blood with CPD has successfully and safely prolonged the time that the patient can be off CRRT
- No adverse events have been reported at our facility on the use of CPD in this method
- Also saves the platelets and plasma in the CRRT circuit blood

MAVO CLINIC

Results and Conclusions

- This procedure did greatly reduce the need to use donor blood to prime CRRT circuits on the patients away from the ICU for less than 8 hours.
- Also been used to save the ECMO circuit blood for reinfusion of pediatric and neonate patients following a successful ECMO wean
- Also this technique was used to save a CRRT circuits blood on a non ECMO pediatric patient



References and Literature

- Maynard, K. (2014). Administration of Blood Components. *AABB Technical Manual* (18th ed., pp. 545-559). Bethesda, MD: AABB.
- Smith, A., Hardison, D., Bridges, B., & Pietsch, J. (2012). Red blood cell transfusion volume and mortality. *Perfusion*, 28(1), 54-60.
- Kochhar, P. K., Ghosh, P., & Kochhar, R. S. (2012). Effect of blood transfusion on subsequent organ transplantation. *Blood Transfusion in Clinical Practice*, 245-254.
- Murphy HJ, Cahill JB, Twombly KE, et al. Implementing a practice change: early initiation of continuous renal replacement therapy during neonatal extracorporeal life support standardizes care and improves short-term outcomes. *J Artif Organs*. 2018 Mar;21(1):76-85.



Acknowledgements

- Thanks to Blake Frazier, BSN, RN, CCRN; Brittany Krom, BSN, RN, CCRN; for working on presentations of this process
- Thank you to the multidisciplinary team that cares of our ECMO patients at Mayo



Questions

- and Thank you

email

- neal.james@mayo.edu

